



cbbL1 Antibody

Product Code	CSB-PA338471XA01FSK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22849
Immunogen	Recombinant Allochromatium vinosum (Chromatium vinosum) cbbL1 protein
Raised In	Rabbit
Species Reactivity	Allochromatium vinosum (Chromatium vinosum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Allochromatium vinosum (Chromatium vinosum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	cbbL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.