



rcsB Antibody

Product Code	CSB-PA361497XA01FUJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05339
Immunogen	Recombinant Klebsiella aerogenes (Enterobacter aerogenes) rcsB protein
Raised In	Rabbit
Species Reactivity	Klebsiella aerogenes (Enterobacter aerogenes)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Klebsiella aerogenes (Enterobacter aerogenes)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rcsB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.