



rbmH Antibody

Product Code	CSB-PA217090XA01LDI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8KI76
Immunogen	Recombinant Lechevalieria aerocolonigenes (Nocardia aerocolonigenes) (Saccharothrix aerocolonigenes) rbmH protein
Raised In	Rabbit
Species Reactivity	Lechevalieria aerocolonigenes (Nocardia aerocolonigenes) (Saccharothrix aerocolonigenes)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Lechevalieria aerocolonigenes (Nocardia aerocolonigenes) (Saccharothrix aerocolonigenes)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rbmH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.