



rplK Antibody

Product Code	CSB-PA320611XA01EZA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10055
Immunogen	Recombinant Proteus vulgaris rplK protein
Raised In	Rabbit
Species Reactivity	Proteus vulgaris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Proteus vulgaris
Lead Time	Made-to-order (14-16 weeks)
Target Names	rplK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.