



rpmA Antibody

Product Code	CSB-PA310172XA01FUL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P94689
Immunogen	Recombinant Prosthecochloris vibrioformis (Chlorobium vibrioforme) rpmA protein
Raised In	Rabbit
Species Reactivity	Prosthecochloris vibrioformis (Chlorobium vibrioforme)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Prosthecochloris vibrioformis (Chlorobium vibrioforme)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rpmA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.