



rpl19e Antibody

Product Code	CSB-PA318460XA01MEZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14024
Immunogen	Recombinant Methanococcus vannielii rpl19e protein
Raised In	Rabbit
Species Reactivity	Methanococcus vannielii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Methanococcus vannielii
Lead Time	Made-to-order (14-16 weeks)
Target Names	rpl19e
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.