



ribH Antibody

Product Code	CSB-PA530305XA01FWR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O68250
Immunogen	Recombinant Sulfurospirillum multivorans (Dehalospirillum multivorans) ribH protein
Raised In	Rabbit
Species Reactivity	Sulfurospirillum multivorans (Dehalospirillum multivorans)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sulfurospirillum multivorans (Dehalospirillum multivorans)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ribH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.