



rpl15 Antibody

Product Code	CSB-PA321191XA01MEZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14032
Immunogen	Recombinant Methanococcus vannielii rpl15 protein
Raised In	Rabbit
Species Reactivity	Methanococcus vannielii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Methanococcus vannielii
Lead Time	Made-to-order (14-16 weeks)
Target Names	rpl15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.