



rnfB Antibody

Product Code	CSB-PA316977XA01RLC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CZ14
Immunogen	Recombinant Rhodobacter capsulatus (Rhodopseudomonas capsulata) rnfB protein
Raised In	Rabbit
Species Reactivity	Rhodobacter capsulatus (Rhodopseudomonas capsulata)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rhodobacter capsulatus (Rhodopseudomonas capsulata)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rnfB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.