



rpl35ae Antibody

Product Code	CSB-PA322623XA01FHZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20299
Immunogen	Recombinant Pyrococcus woesei rpl35ae protein
Raised In	Rabbit
Species Reactivity	Pyrococcus woesei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pyrococcus woesei
Lead Time	Made-to-order (14-16 weeks)
Target Names	rpl35ae
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.