



rpmH Antibody

Product Code	CSB-PA358015XA01FFZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A162
Immunogen	Recombinant Pseudomonas putida (Arthrobacter siderocapsulatus) rpmH protein
Raised In	Rabbit
Species Reactivity	Pseudomonas putida (Arthrobacter siderocapsulatus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudomonas putida (Arthrobacter siderocapsulatus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rpmH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.