



rnhA Antibody

Product Code	CSB-PA359051XA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A7Y4
Immunogen	Recombinant Escherichia coli (strain K12) rnhA protein
Raised In	Rabbit
Species Reactivity	Escherichia coli (strain K12)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli (strain K12)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rnhA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.