



RRF Antibody

Product Code	CSB-PA307806XA01FKI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P82231
Immunogen	Recombinant Spinacia oleracea (Spinach) RRF protein
Raised In	Rabbit
Species Reactivity	Spinacia oleracea (Spinach)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Spinacia oleracea (Spinach)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	RRF