



rpoC2 Antibody

Product Code	CSB-PA295989XA01NEF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14564
Immunogen	Recombinant Nostoc commune rpoC2 protein
Raised In	Rabbit
Species Reactivity	Nostoc commune
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Nostoc commune
Lead Time	Made-to-order (14-16 weeks)
Target Names	rpoC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.