



RRF Antibody

Product Code	CSB-PA334308XA01FWV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37706
Immunogen	Recombinant Daucus carota (Wild carrot) RRF protein
Raised In	Rabbit
Species Reactivity	Daucus carota (Wild carrot)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Daucus carota (Wild carrot)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	RRF