



rpsL Antibody

Product Code	CSB-PA362830XA01GFM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09901
Immunogen	Recombinant Geobacillus stearothermophilus (Bacillus stearothermophilus) rpsL protein
Raised In	Rabbit
Species Reactivity	Geobacillus stearothermophilus (Bacillus stearothermophilus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Geobacillus stearothermophilus (Bacillus stearothermophilus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rpsL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.