



RNH1 Antibody

Product Code	CSB-PA248172XA01DXW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07762
Immunogen	Recombinant Crithidia fasciculata RNH1 protein
Raised In	Rabbit
Species Reactivity	Crithidia fasciculata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Crithidia fasciculata
Lead Time	Made-to-order (14-16 weeks)
Target Names	RNH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.