



rpoB Antibody

Product Code	CSB-PA888654XA01FWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9KW14
Immunogen	Recombinant Shewanella violacea rpoB protein
Raised In	Rabbit
Species Reactivity	Shewanella violacea
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Shewanella violacea
Lead Time	Made-to-order (14-16 weeks)
Gene Names	rpoB