



rtpR Antibody

Product Code	CSB-PA682405XA01LAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q59490
Immunogen	Recombinant Lactobacillus leichmannii rtpR protein
Raised In	Rabbit
Species Reactivity	Lactobacillus leichmannii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Lactobacillus leichmannii
Lead Time	Made-to-order (14-16 weeks)
Target Names	rtpR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.