



rseP Antibody

Product Code	CSB-PA882770XA01PII
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9S342
Immunogen	Recombinant Photorhabdus luminescens (Xenorhabdus luminescens) rseP protein
Raised In	Rabbit
Species Reactivity	Photorhabdus luminescens (Xenorhabdus luminescens)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Photorhabdus luminescens (Xenorhabdus luminescens)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rseP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.