



sbmC Antibody

Product Code	CSB-PA509911XA01DTP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D2TPS8
Immunogen	Recombinant Citrobacter rodentium (strain ICC168) (Citrobacter freundii biotype 4280) sbmC protein
Raised In	Rabbit
Species Reactivity	Citrobacter rodentium (strain ICC168) (Citrobacter freundii biotype 4280)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Citrobacter rodentium (strain ICC168) (Citrobacter freundii biotype 4280)
Lead Time	Made-to-order (14-16 weeks)
Target Names	sbmC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.