



PFE1050w Antibody

Product Code	CSB-PA001474XA01EWP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50250
Immunogen	Recombinant Plasmodium falciparum (isolate 3D7) PFE1050w protein
Raised In	Rabbit
Species Reactivity	Plasmodium falciparum (isolate 3D7)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Plasmodium falciparum (isolate 3D7)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	PFE1050w