



Isc Antibody

Product Code	CSB-PA525146XA01FWZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O52408
Immunogen	Recombinant Pseudomonas savastanoi pv. glycinea (Pseudomonas syringae pv. glycinea) Isc protein
Raised In	Rabbit
Species Reactivity	Pseudomonas savastanoi pv. glycinea (Pseudomonas syringae pv. glycinea)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudomonas savastanoi pv. glycinea (Pseudomonas syringae pv. glycinea)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Isc
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.