



Toxin BmKT Antibody

Product Code	CSB-PA839222XA01MRE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q95P69
Immunogen	Recombinant Mesobuthus martensii (Manchurian scorpion) (Buthus martensii) Toxin BmKT protein
Raised In	Rabbit
Species Reactivity	Mesobuthus martensii (Manchurian scorpion) (Buthus martensii)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mesobuthus martensii (Manchurian scorpion) (Buthus martensii)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.