



Uncharacterized 9.2 kDa Antibody

Product Code	CSB-PA325647XA01FQX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19195
Immunogen	Recombinant Escherichia phage Bf23 (Enterobacteria phage BF23) Uncharacterized 9.2 kDa protein
Raised In	Rabbit
Species Reactivity	Escherichia phage Bf23 (Enterobacteria phage BF23)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage Bf23 (Enterobacteria phage BF23)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.