



orf2 Antibody

Product Code	CSB-PA313878XA01FMW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C2G9
Immunogen	Recombinant Streptococcus pneumoniae orf2 protein
Raised In	Rabbit
Species Reactivity	Streptococcus pneumoniae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptococcus pneumoniae
Lead Time	Made-to-order (14-16 weeks)
Gene Names	orf2