



# HI\_0278 Antibody

<b>Product Code</b>	CSB-PA337599XA01HTA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P43976
<b>Immunogen</b>	Recombinant Haemophilus influenzae (strain ATCC 51907 / DSM 11121 / KW20 / Rd) HI_0278 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Haemophilus influenzae (strain ATCC 51907 / DSM 11121 / KW20 / Rd)
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Haemophilus influenzae (strain ATCC 51907 / DSM 11121 / KW20 / Rd)
<b>Lead Time</b>	Made-to-order (14-16 weeks)
<b>Target Names</b>	HI_0278
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.