



SYNPCC7002_A0425 Antibody

Product Code	CSB-PA324847XA01FQA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19737
Immunogen	Recombinant Synechococcus sp. (strain ATCC 27264 / PCC 7002 / PR-6) (Agmenellum quadruplicatum) SYNPCC7002_A0425 protein
Raised In	Rabbit
Species Reactivity	Synechococcus sp. (strain ATCC 27264 / PCC 7002 / PR-6) (Agmenellum quadruplicatum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Synechococcus sp. (strain ATCC 27264 / PCC 7002 / PR-6) (Agmenellum quadruplicatum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	SYNPCC7002_A0425
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.