



xpaL2 Antibody

Product Code	CSB-PA859916XA01BQT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99165
Immunogen	Recombinant Bacillus licheniformis xpaL2 protein
Raised In	Rabbit
Species Reactivity	Bacillus licheniformis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus licheniformis
Lead Time	Made-to-order (14-16 weeks)
Target Names	xpaL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.