



Uncharacterized protein UU427 homolog Antibody

Product Code	CSB-PA690905XA01UAP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q56564
Immunogen	Recombinant Ureaplasma urealyticum (Ureaplasma urealyticum biotype 2) Uncharacterized protein UU427 homolog protein
Raised In	Rabbit
Species Reactivity	Ureaplasma urealyticum (Ureaplasma urealyticum biotype 2)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Ureaplasma urealyticum (Ureaplasma urealyticum biotype 2)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.