



Uncharacterized protein ORF43.1 in regA 3'region Antibody

Product Code	CSB-PA273448XA01EDR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18244
Immunogen	Recombinant Enterobacteria phage RB51 (Bacteriophage RB51) Uncharacterized protein ORF43.1 in regA 3'region protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage RB51 (Bacteriophage RB51)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage RB51 (Bacteriophage RB51)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.