



yidC Antibody

Product Code	CSB-PA328352XA01EYY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22833
Immunogen	Recombinant Proteus mirabilis yidC protein
Raised In	Rabbit
Species Reactivity	Proteus mirabilis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Proteus mirabilis
Lead Time	Made-to-order (14-16 weeks)
Target Names	yidC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.