



xpaL1 Antibody

Product Code	CSB-PA857849XA01BQT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99164
Immunogen	Recombinant Bacillus licheniformis xpaL1 protein
Raised In	Rabbit
Species Reactivity	Bacillus licheniformis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus licheniformis
Lead Time	Made-to-order (14-16 weeks)
Gene Names	xpaL1