



yebY Antibody

Product Code	CSB-PA353120XA01SZB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P64507
Immunogen	Recombinant Shigella flexneri yebY protein
Raised In	Rabbit
Species Reactivity	Shigella flexneri
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Shigella flexneri
Lead Time	Made-to-order (14-16 weeks)
Target Names	yebY
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.