



# yfeD Antibody

Product Code	CSB-PA690921XA01YAS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q56955
Immunogen	Recombinant Yersinia pestis yfeD protein
Raised In	Rabbit
Species Reactivity	Yersinia pestis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Yersinia pestis
Lead Time	Made-to-order (14-16 weeks)
Gene Names	yfeD