



Uncharacterized family 31 glucosidase ORF2 Antibody

Product Code	CSB-PA312454XA01EOS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01336
Immunogen	Recombinant Escherichia vulneris Uncharacterized family 31 glucosidase ORF2 protein
Raised In	Rabbit
Species Reactivity	Escherichia vulneris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia vulneris (Pseudoescherichia vulneris)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.