



Uncharacterized 16.9 kDa protein in acs 3'region Antibody

Product Code	CSB-PA104107XA01FRD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27096
Immunogen	Recombinant Methanosaeta concilii (Methanotherix soehngeni) Uncharacterized 16.9 kDa protein in acs 3'region protein
Raised In	Rabbit
Species Reactivity	Methanosaeta concilii (Methanotherix soehngeni)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Methanotherix soehngeni (Methanosaeta concilii)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.