



# ybbY Antibody

|                     |                                                                                  |
|---------------------|----------------------------------------------------------------------------------|
| Product Code        | CSB-PA300585XA01ENV                                                              |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | P77328                                                                           |
| Immunogen           | Recombinant Escherichia coli (strain K12) ybbY protein                           |
| Raised In           | Rabbit                                                                           |
| Species Reactivity  | Escherichia coli (strain K12)                                                    |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                     |
| Form                | Liquid                                                                           |
| Conjugate           | Non-conjugated                                                                   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified                                                        |
| Isotype             | IgG                                                                              |
| Clonality           | Polyclonal                                                                       |
| Product Type        | Polyclonal Antibody                                                              |
| Immunogen Species   | Escherichia coli (strain K12)                                                    |
| Lead Time           | Made-to-order (14-16 weeks)                                                      |
| Gene Names          | ybbY                                                                             |