



Uncharacterized protein in rusA 5'region Antibody

Product Code	CSB-PA656412XA01ECB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q37870
Immunogen	Recombinant Enterobacteria phage 82 (Bacteriophage 82) Uncharacterized protein in rusA 5'region protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage 82 (Bacteriophage 82)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage 82 (Bacteriophage 82)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.