



Uncharacterized protein Gp1 Antibody

Product Code	CSB-PA323348XA01EDB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15234
Immunogen	Recombinant Enterobacteria phage N4 (Bacteriophage N4) Uncharacterized protein Gp1 protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage N4 (Bacteriophage N4)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage N4 (Bacteriophage N4)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.