



Uncharacterized protein in acs 3'region Antibody

Product Code	CSB-PA969993XA01FRD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27097
Immunogen	Recombinant Methanosaeta concilii (Methanotherix soehngeni) Uncharacterized protein in acs 3'region protein
Raised In	Rabbit
Species Reactivity	Methanosaeta concilii (Methanotherix soehngeni)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Methanotherix soehngeni (Methanosaeta concilii)
Lead Time	Made-to-order (14-16 weeks)