



Uncharacterized protein in nthA 5'region Antibody

Product Code	CSB-PA247080XA01RLJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20771
Immunogen	Recombinant Rhodococcus erythropolis (Arthrobacter picolinophilus) Uncharacterized protein in nthA 5'region protein
Raised In	Rabbit
Species Reactivity	Rhodococcus erythropolis (Arthrobacter picolinophilus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rhodococcus erythropolis (Arthrobacter picolinophilus)
Lead Time	Made-to-order (14-16 weeks)