



UPF0194 membrane protein in asrC 5'region Antibody

Product Code	CSB-PA867902XA01FRQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9L9D4
Immunogen	Recombinant Acidithiobacillus ferrooxidans (Thiobacillus ferrooxidans) UPF0194 membrane protein in asrC 5'region protein
Raised In	Rabbit
Species Reactivity	Acidithiobacillus ferrooxidans (Thiobacillus ferrooxidans)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Acidithiobacillus ferrooxidans (Thiobacillus ferrooxidans)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.