



Extracellular serine proteinase Antibody

Product Code	CSB-PA303605XA01TDV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80146
Immunogen	Recombinant Thermus sp. (strain Rt41A) Extracellular serine proteinase protein
Raised In	Rabbit
Species Reactivity	Thermus sp. (strain Rt41A)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Thermus sp. (strain Rt41A)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.