



SHBG Antibody

Product Code	CSB-PA021252XA01RB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15196
Immunogen	Recombinant Oryctolagus cuniculus (Rabbit) SHBG protein
Raised In	Rabbit
Species Reactivity	Oryctolagus cuniculus (Rabbit)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Lead Time	Made-to-order (14-16 weeks)
Target Names	SHBG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.