



sotB Antibody

Product Code	CSB-PA879661XA01FRO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9S3J9
Immunogen	Recombinant Dickeya chrysanthemi (Pectobacterium chrysanthemi) (Erwinia chrysanthemi) sotB protein
Raised In	Rabbit
Species Reactivity	Dickeya chrysanthemi (Pectobacterium chrysanthemi) (Erwinia chrysanthemi)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Dickeya chrysanthemi (Pectobacterium chrysanthemi) (Erwinia chrysanthemi)
Lead Time	Made-to-order (14-16 weeks)
Target Names	sotB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.