



sdhA Antibody

Product Code	CSB-PA339321XA01FWX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33073
Immunogen	Recombinant Peptoniphilus asaccharolyticus (Peptostreptococcus asaccharolyticus) sdhA protein
Raised In	Rabbit
Species Reactivity	Peptoniphilus asaccharolyticus (Peptostreptococcus asaccharolyticus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Peptoniphilus asaccharolyticus (Peptostreptococcus asaccharolyticus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	sdhA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.