



shc Antibody

Product Code	CSB-PA335814XA01FXF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33247
Immunogen	Recombinant Alicyclobacillus acidocaldarius subsp. acidocaldarius (strain ATCC 27009 / DSM 446 / JCM 5260 / NBRC 15652 / NCIMB 11725 / NRRL B-14509 / 104-1A) (Bacillus acidocaldarius) shc protein
Raised In	Rabbit
Species Reactivity	Alicyclobacillus acidocaldarius subsp. acidocaldarius (strain ATCC 27009 / DSM 446 / JCM 5260 / NBRC 15652 / NCIMB 11725 / NRRL B-14509 / 104-1A) (Bacillus acidocaldarius)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Alicyclobacillus acidocaldarius subsp. acidocaldarius (strain ATCC 27009 / DSM 446 / JCM 5260 / NBRC 15652 / NCIMB 11725 / NRRL B-14509 / 104-1A) (Bacillus acidocaldarius)
Lead Time	Made-to-order (14-16 weeks)
Target Names	shc
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.