



GH1 Antibody

Product Code	CSB-PA816827XA01AYX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8HYE5
Immunogen	Recombinant Ailuropoda melanoleuca (Giant panda) GH1 protein
Raised In	Rabbit
Species Reactivity	Ailuropoda melanoleuca (Giant panda)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Ailuropoda melanoleuca (Giant panda)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.