



Flavonoid 4'-O-methyltransferase Antibody

Product Code	CSB-PA547588XA01GGV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C6TAY1
Immunogen	Recombinant Glycine max (Soybean) (Glycine hispida) Flavonoid 4'-O-methyltransferase protein
Raised In	Rabbit
Species Reactivity	Glycine max (Soybean) (Glycine hispida)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Glycine max (Soybean) (Glycine hispida)
Lead Time	Made-to-order (14-16 weeks)